Best. Nr. CA 015
Inhalt: 1 x 100 mL Farbreagenz
 1 x 100 mL Pufferlösung

Zusätzlich erforderlich:
 Calcium-Magnesium-Standard Best. Nr. CA ST
 Diaglobal-Rundküvetten Best. Nr. LH 075

Methode
 o-CPC-Methode¹⁾

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Reagent for quantitative In-vitro-determination of Calcium in plasma / urine of animals

Order No. CA 015
Contents 1 x 100 mL Colour reagent
 1 x 100 mL Buffer solution
Additionally needed:
 Calcium-Magnesium-Standard Order no. CA ST
 Diaglobal round cuvettes Order no. LH 075

Method
 o-CPC-method¹⁾

Sample material
 Heparinised plasma or urine
 Do not use EDTA-plasma
 Stability in plasma: at +2°C to +25°C: 10 days
 Stability in urine: at +2°C to +8°C: 7 days
 Sedimentation of calciumsalts can be prevented by adding of 10 mL concentrated hydrochloric acid.

Reagent
 Contents / concentrations:
 1. Colour reagent
 o-Cresolphthalein-complexone 0.16 mmol/L;
 Hydroxyquinoline 7 mmol/L
 2. Buffer solution
 Glycine 300 mmol/L; Sodium azide < 0.1%;
 Triton X-100 < 1%

Safety information
 The reagents contain sodium azide (< 0.1 %) and Triton X-100 (< 1%). Do not swallow and avoid contact with skin and mucous membranes.
 If desired a safety data sheet will be provided.²⁾

Storage and shelf life
 The reagents can be kept at a temperature between +15°C and +25°C until the expiry date indicated on the packaging.

Measurement conditions
 Measurement device: VET Photometer Diaglobal
 Meas. wavelength: 520 nm
 Temperature: Room temperature

Measurement range
 0.10 - 5.00 mmol/L (0.40 - 20.0 mg/dL)
 Should values exceed this range, dilute sample 1+1 with physiological saline solution. Multiply the result by 2.
 Dilute urine 1+2 (result x 3) or 1+4 (result x 5)

Notes
 1. The calcium assay is very sensitive. We recommend the application of disposables.
 2. For the "on-site" preparation of plasma we recommend the Diaglobal - Mini centrifuge (Order No. DZ 002)

Working instructions

Pipette into round cuvettes:			
	Blank	Standard	Analysis
Colour reagent	1500 µL	1500 µL	1500 µL
Buffer solution	1500 µL	1500 µL	1500 µL
Standard	--	50 µL	--
Sample	--	--	50 µL

Mix well and measure immediately or within 1h standard and analysis against blank.

- Select the <CA> test
- Insert cuvette with reagent into photometer (blank value)
- Remove blank cuvette after signal ton
- Insert cuvette with standard into photometer (standard-absorption is shown on the display)
- Remove cuvette
- Insert cuvette with sample (Analysis) into photometer
- Read the result

Tip
 It is possible to prepare several vials for measurement at the same time by pre-pipetting 1500 µL of staining medium and 1500 µL of buffer solution into the appropriate number of Diaglobal round cuvettes. Shelf life of the pipetted reagents in the vial: 5 hours

Calculation
 Concentration of calcium in plasma/urine:
 $c \text{ [mmol/L]} = 2.50 \times \text{Abs. analysis} / \text{Abs. standard}$
 $c \text{ [mg/dL]} = 10 \times \text{Abs. analysis} / \text{Abs. standard}$
 Conversion in mg/dL: mmol/L x 4.008 = mg/dL

Calculation of release of calcium per day:
 Concentration of calcium [mmol/L] x 24h-Vol. of urine [L] = mmol calcium/24 h

Quality assurance
 For quality assurance we recommend our calcium - /magnesium - control solution CA QS.

Reference values for cattle


	mmol/L	mg/dL
Plasma	2.00 - 3.00	8.00 - 12.0

Summary³⁾
 Approximately 99% of the calcium in mammalian organisms is localised in the bone structures. The remainder of the calcium is distributed around the other tissues and extracellular fluids, where it is important for, among other things, the functioning of the muscles and the nervous system, lactation and digestion.
 The calcium level in the plasma is key to controlling calcium metabolism. Fifty-five percent of the calcium in the plasma is ionised calcium (Ca²⁺), 40% is bound to protein and 5% to organic acids. The total calcium is determined using the Diaglobal CA 015 test.

Significance
 Determining the level of calcium plays a key role in the diagnosis, prophylaxis and treatment of milk fever in cows - a disease which occurs shortly after calving and is characterised by low calcium levels.
 Even calcium levels between 1.75 and 2.00 mmol/L lead to typical symptoms (subclinical milk fever) which require an intervention.

At levels <1.75 mmol/L, acute hypocalcaemia occurs (clinical milk fever), which leads to the inability of the cow to stand and is associated with symptoms such as apathy, low temperature, low rumen activity and reduction in the milk output.
 High-yield cows which require increased levels of calcium due to their high milk production are at particular risk.
 The Diaglobal CA 015 test was designed for "on-site" use. The plasma needed for the determination can be extracted using a mobile mini centrifuge.

Measurement principle⁴⁾
 In alkaline buffer solution, calcium forms with o-cresolphthalein-complexone (o-CPC) a violet dye.

Adding 8-hydroxyquinoline to the reagent avoids interferences by magnesium-ions and iron-ions.
 The colour intensity of Calcium-o-CPC is directly proportional to the concentration of calcium in serum. Calculation is due by standard.

Performance parameters
Specificity / interferences
 No significant impact caused by haemoglobin (< 2.0 g/L) and bilirubin (< 20 mg/dL). Severe lipidaemia interferes. Complexing anticoagulants (EDTA, citrate and oxalate) may not be used.
 Weak complexing agents such as gluconic acid and L-lactic acid, which are sometimes used as calcium salts in animal feed, do not interfere.

Inaccuracy
 The reproducibility was checked using control samples, human serums and cow plasmas. Typical VK values are shown below.

Human serum

In series [n = 20]	Average [mmol/L]	Standard deviation [mmol/L]	VK [%]
Sample 1	2.02	0.04	2.0
Sample 2	2.69	0.02	0.9
Sample 3	3.57	0.06	1.7

From day to day [n = 20]	Average [mmol/L]	Standard-deviation [mmol/L]	VK [%]
Sample 1	1.77	0.05	2.7
Sample 2	2.41	0.06	2.3
Sample 3	3.66	0.08	2.1

Cattle plasma

In series [n = 20]	Average [mmol/L]	Standard-deviation [mmol/L]	VK [%]
Sample 1	1.66	0.02	1.4
Sample 2	2.31	0.04	1.8

Analytic sensitiveness
 Lower detection limit: 0.10 mmol/L (0.40 mg/dL)

Comparison of methods
 A comparison of the Diaglobal test CA 015 (y) and the original methode (x) (Spectral photometer, measurement at 578nm) resulted in the following correlation according to the linear regression:
 $y = 0.972x - 0.015$
 $r = 0.998$

Sample material: heparinised cattle plasma
 n = 20, concentration range: 0.8 - 3.5 mmol/L

Information on disposal
 Waste code number 180106:
 Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system.
 Dispose of in accordance with official regulations.
 Non-contaminated and completely empty packaging can be recycled.

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