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Order No. CHO 142
Content: 40 tests

Method

Enzymatic colorimetric test, CHOD-PAP method¹⁾
The determination can directly employ both whole blood and serum/plasma. The blood is immediately and completely hemolysed by the reagent. The result obtained from blood samples refers to the plasma concentrations.

Sample material

Capillary blood or EDTA venous blood
Serum, heparinized or EDTA plasma
Pipette capillary blood immediately in single test cuvettes.
Stability of the haemolyzed sample in the buffer solution:
at +2 to +8°C: 8 hours
at +15 to +30°C: 4 hours

Reagent

Contents / concentrations:
1. Starter reagent (caps in PE-bottle)
Cholesteroloxidase (CHOD) from Brevi bacterium >350 U/L,
Peroxidase (POD) >4 kU/L, 4-Aminophenazone 0.20 mmol/L
2. Buffer solution (pre-portioned in round cuvettes)
Lipoprotein lipase / Cholesterol esterase from microorg. >800 U/L,
4-Chlorophenol 9 mmol/L, Sodium azide <0.1 %, Triton X-100 < 1%, PIPES buffer 75 mmol/L, pH 7.6

Safety information

The buffer solution (round cuvette) contains Sodium azide (<0.1 %) and Triton X-100. Do not swallow and avoid contact with skin and mucous membranes. If desired a safety data sheet will be provided.²⁾

Storage and shelf life

Reagents can be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging. Please take the screw caps out of the container just before the analysis and close the container immediately.

Measurement conditions

Measurement devices: Diaglobal Photometer

Meas. wavelength: 520nm
Temperature: Room temperature

Measurement range

20 - 1000 mg/dL (0.52 - 25.9 mmol/L)

Working instructions

The measurement can be performed as a single or serial measurement [with a balancing of the A(0) counts].

A. Single measurement

Pipette in round cuvette:	
	Analysis
Sample	10 µL
Mix thoroughly.	

- Select the <CHO> test
- Insert analysis cuvette (blank value)
- Screw the cap from PE-bottle onto the cuvette, dissolve the starter reagent by inverting several times
- Press [ON/ENTER]
- Insert analysis cuvette again
- Wait for result

B. Measurement of series (up to 20 samples)

Pipette in round cuvettes:	
	Analysis
Sample	10 µL
Mix thoroughly.	

- Select the <CHO> test
- Insert the analysis cuvettes one after another (blank values)
- Screw the caps from PE-bottle onto the cuvettes, dissolve the starter reagent by inverting several times
- Press [ON/ENTER]
- Insert the first analysis cuvette **immediately** again
- Wait for result
- Insert the other analysis cuvettes one after another in the same order as of the blank value measurement
- Results of the respective analysis cuvette can be read immediately

Quality assurance

For quality assurance we recommend universal control sera from company Roche, www.roche.de:
PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
Order-No.: 05 947 626 190 / 05 947 774 190
Ref.: Roche / Hitachi analyzers,
Method: CHOD-PAP stand. ID/MS

Clinical interpretation

The following values are recommended to assess the risk factor of cholesteremia:³⁾

	mg/dL	mmol/L
Low risk	> 200	> 5.2
High risk	> 260	> 6.7

Summary

Cholesterol pertains to the essential elements of cell membranes. Furthermore, it is the initial substance for the synthesis of steroid hormones and bile acids. Circa 25 % derive from nutrition. The rest is synthesised in the organism. Cholesterol is carried in the plasma as a lipoprotein complex due to its inability to solubilise in water.

Indications / diagnostic significance:¹⁾


- Early diagnosis of the risk of atherosclerosis
- Control of lipid lowering therapies

The determination of cholesterol as a whole is counted among the basis programme of cardiovascular preventive medical checkups. According to the recommendations of the European Atherosclerosis Society, a cholesterol as a whole count of 200 mg/dL is classified as marginal value whose exceeding initially means a moderately, then intensely arisen risk of atherosclerosis.³⁾ This risk increases because of factors like smoking and hypertension. The determination of the cholesterol proportions, which are associated with the lipoprotein fractions HDL and LDL, allows further risk evaluation.

In order to determine cholesterol as a whole the CHOD-PAP method is applied without exception nowadays. This method forms the basis of the Diaglobal test. The non-enzymatical Liebermann-Burchard reaction in turn has become irrelevant for routine diagnostics.

Measurement principle

The cholesterol esters are decomposed enzymatically. After this the free cholesterol becomes oxidised by means of the cholesteroloxidase (CHOD) and the generated H₂O₂ is converted into a dye through a PAP reaction.⁴⁾

The concentration of the quinonimine dye is a measure for the cholesterol concentration in the sample and is measured photometrically. The end point of the reaction is identified automatically by the photometer.

Performance parameters

Specificity / interferences⁵⁾

Too high values because of haemoglobin (> 2.0 g/L), no interference due to lipaemia, ascorbic acid in physiological

concentrations (< 50 mg/dL) and bilirubin (< 15 mg/dL).
Pharmaceutical interferences: low values due to methyldopa (> 50 mg/L) and novaminsulfon (>150 mg/L) in concentrations above the therapeutic field.

Inaccuracy

The reproducibility was checked using human and control samples.

In series [n = 20]	Average [mg/dL]	Standard deviation [mg/dL]	VK [%]
EDTA blood	188	5	2.5
Serum 1	201	3	1.5
Serum 2	402	5	1.3
From day to day [n = 20]	Average [mg/dL]	Standard deviation [mg/dL]	VK [%]
Serum 1	199	3	1.6
Serum 2	398	6	1.5

Analytic sensitiveness

Lower detection limit: 20 mg/dL (0.52 mmol/L)

Comparison of methods

A comparison of the Diaglobal test CHO 142 (y) and a commercially available test (x) based on the CHOD-PAP method as well as a comparison between values which were determined with blood (y) and plasma (x) of the respectively same proband resulted in the following correlation data according to the Passing/Bablok⁶⁾ process:

a) plasma
 $y = 0,990x + 3,5$
 $r = 0,997$
 $n = 60$

b) EDTA blood / plasma
 $y = 1,000x + 1,0$
 $r = 0,988$
 $n = 85$

Concentration range: 94 - 501 mg/dL

Information on disposal

Waste code number 180106:
Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system. Dispose of in accordance with official regulations. Non-contaminated and completely empty packaging can be recycled.

Bibliography

1. Thomas L. Labor und Diagnose. 4th edition. Marburg: Die Medizinische Verlagsgesellschaft 1995: 200
2. <http://www.diaglobal.de/de/service/downloads/index.html>
3. Study Group, European Atherosclerosis Society. Eur Heart J. 1988; The recognition and management of hyperlipidemia in adults.
4. Allain CC. Clin Chem 1974; 20:470
5. Sonntag O. Arzneimittelinterferenzen. Stuttgart, New York: Thieme Verlag, 1985:24
6. Passing H, Bablok W. A new biometric procedure for testing the equality of measurements from two different analytical methods. J Clin Chem Clin Biochem. 1983; 21:709-720