**Best. Nr.** GLY 742  
**Inhalt:** 40 Tests

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
- <http://www.diaglobal.de/de/service/downloads/index.html>

**Order No.** GLY 742  
**Contents:** 40 tests

**Method**  
Enzymatic colorimetric test, GPO-PAP method

**Reagents**  
Contents / concentrations of the ready-to-use solution:

1. Starter reagent (screw caps)
  - L-glycerol-3-phosphate-oxidase (GPO) from microorganisms > 3.5 kU/L, Glycerokinase (GK) from bacillus stearothermophilus > 0.9 kU/L, Peroxidase (POD) >3.5 kU/L, ATP 2.4 mmol/L, 4-Aminophenazone 0.15 mmol/L
2. Buffer solution (pre-portioned in round cuvettes)
  - 2,4-Dichlorophenol 4 mmol/L, Sodium azide < 0.1 %, Triton X-100 < 1%, PIPES-buffer 50 mmol/L, pH 7.5

**Safety information**  
The buffer solution (round cuvette) contains sodium azide (< 0.1 %) and Triton X-100. Do not swallow and avoid contact with skin and mucous membranes. If desired a safety data sheet will be provided.<sup>1)</sup>

**Storage and shelf life**  
The test reagents can be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging. Please take the screw caps out of the container just before the analysis.

**Measurement conditions**  
Measurement device: Biodiesel Photometer Diaglobal

Meas. wavelength: 520nm

Temperature: Room temperature

The algorithm to compute the analysis result is coded in the above-named photometer.

### Working instructions for aqueous solutions

**Sample material**  
Aqueous glycerol solutions

**A. Compounds of glycerol**  
Dilute sample 1 : 1000 in this way: Weigh exactly 1.26 g from the glycerol compound, dilute with dest. water and fill up to 1 L in volumetric flask. Take 10 µL from this solution for measurement.  
The measurement can be performed as a single as well as a serial measurement (with a balancing of the A (0) values = blank values).

Pipette into round cuvette:	
	Analysis
Sample	10 µL
Mix thoroughly	

- Select <GLY aqu.>
- Insert analysis cuvette (blank value)
- Screw the cap from PE-bottle onto the cuvette, dissolve the starter reagent by inverting several times
- Press [ON/ENTER]
- Insert cuvette into the photometer immediately
- After the end of the reaction the measurement result of the sample is shown in mg/dL and in % resp.

**Calculation**  
Result on display = GLY [mg/dL]  
= GLY [%]

**Measurement range**  
12.6 - 200 mg/dL resp. 10 - 110 %

**B. Aqueous glycerol solutions**  
Fill up 1 mL solution to 1 L (1:1000) with dest. water. Take 10 µL from this solution for measurement. Procedure of measurement like stated above.  
If the glycerol concentration is expected to be less than 10%, please dilute as follows:  
Dilute sample 1: 100, result / 10  
Dilute sample 1: 10, result / 100

**Note**  
For the calculation of the glycerol concentration in percent v/v, the density of glycerol 1.26 g/cm<sup>3</sup> has been used.

### Working instructions for biodiesel

**Sample material**  
Biodiesel

Additionally required:  
**LH 051** Reaction tubes with polystyrene pearls

**Working instructions**  
Indissoluble fatty acid methyl ester which may create a turbidity in buffer solution will be eliminated through the adsorption on polystyrene pearls.

- Pipette **1000 µL** dest. water in reaction tube with polystyrene pearls
- Pipette **20 µL** biodiesel with micropipettor on polystyrene pearls and wash out well
- Mix thoroughly by shaking it well
- Wait 1 minute
- Pipette **500 µL** aqueous solution into cuvette GLY 742 with buffer solution
- Select <GLY org.>
- Insert analysis cuvette (blank value)
- Screw the screw cap from PE-bottle onto the cuvette, dissolve the starter reagent by inverting several times
- Press [ON/ENTER]
- Insert cuvette into the photometer immediately
- At the end of the reaction time (2 - 3 min) the measurement result will be displayed in Ma%

**Calculation:**  
Concentration of glycerol:  
c (mg/dL) = absorbance x factor  
c (Ma %) = c (mg/dL) / density (biodiesel) / 1000  
This value is shown on the display.  
Density of biodiesel: 0,8776 g/cm<sup>3</sup> has been used

**Measurement range**  
0.001 - 0.250 Ma%

**Dilute when exceeding the measurement range**  
In case of exceeding measurement range, on display will shown > 0.25 Ma%. In that case dilute the rest aqueous solution of reaction tube with dest. water (1:20) and repeat the measurement with 500 µL of this dilution.  
Result x 20

**Note**  
If the glycerol concentration required only for correction of triglycerid value please perform the measurement with reagent GLY 742 in mode <TRI>. Result in g/dL TRI

**Information on disposal**  
Waste code number 180106:  
Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system. Dispose of in accordance with official regulations. Non-contaminated and completely empty packaging can be recycled. Non-contaminated and completely empty packaging can be recycled.

**Bibliography**  
1. <http://www.diaglobal.de/de/service/downloads/index.html>