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Article No. HB 142
Name Haemoglobin Kit (HiCN)
Contents 40 round cuvettes
 Ready-to-use reagent (2.5 mL)

Short description
 Kit for the photometric determination of the concentration of total haemoglobin in the blood according to the cyanmethaemoglobin method (HiCN)

Manual, non-automated working method for the quantitative determination of the hemoglobin concentration with a Diaglobal photometer
 In-vitro diagnostic device for single use

Sample material
 Capillary blood or venous blood

Intended use
 Detection and follow-up of anemia and polyglobulia as well as for diagnosis and monitoring of people at risk from iron deficiency

Suitable for examinations in the laboratory and near-patient testing
 Not intended for self-testing

User
 Suitable for near-patient diagnostics
 Ambulances, hospital emergency rooms, registered doctors
 Gynecologists in the context of prenatal care who carry out haemoglobin measurements to clarify suspected anemia in their own practices
 Suitable for monitoring therapeutic measures such as prescribing iron tablets. The measurement takes place in the presence of the patient, so that a treatment plan can be initiated immediately.

Measuring principle
 The cyanmethaemoglobin method, which is the international reference method according to the recommendations of the CLSI
 By lysing of the erythrocyte membrane, haemoglobin is transferred into the reaction solution, oxidized to methaemoglobin with potassium hexacyanoferrate (III) and converted with cyanide into cyanmethaemoglobin, which has a broad absorption maximum at 540nm.

The absorption of the cyanmethaemoglobin is proportional to the haemoglobin concentration in the sample and is measured photometrically.

Reagent
 Ready-to-use aqueous reagent (2.5 mL) in round cuvettes
 Concentration of the active ingredients:
 Potassium hexacyanoferrate (III) 0.6 mmol/L
 Potassium cyanide 0,7 mmol/L
 Sodium hydrogen carbonate 18 mmol/L
 pH value 8.0

Storage and Shelf Life
 Store at 15 to 25 °C
 Protect from strong light
 Do not freeze
 Best before the expiry date stated on the packaging

Analytical measuring range
 1.0 - 25 g/dL (0.62 - 15.5 mmol/L)

Measuring devices and measuring conditions

Measuring device	Diaglobal photometer
Measuring wavelength	546 nm
Measurement against	Reagent blank value
Measuring temperature	Room temperature

Calibrators
 Diaglobal photometers are factory calibrated; calibration by the user is not possible.

Collection and treatment of the sample
 The extraction is carried out with a ring marked capillary (10 µL)
 Do not use end-to-end capillaries

Use capillary blood immediately for determination. Avoid pressing too hard with the fingertip, otherwise the blood to be extracted will be diluted by tissue fluid.

Whole blood, and venous blood (EDTA or heparin as an anticoagulant) can be stored for up to 24 hours at or below 25°C.
 Sample preparation is not required.

Instructions

Pipette HB 142 into round cuvette	
	Analysis
Blood	10 µL
Eject blood with a micropipette, rinse the capillary by pulling it up and ejecting the sample several times	
Close the round cuvette, mix well, measure after at least 3 min and within 30 min	
Select test <HB>	
Set the zero point with the unprocessed round cuvette, remove the round cuvette after the signal tone	
Place the round cuvette with the blood sample in the photometer	
Read the measured value	

Calculation
 The haemoglobin concentration is calculated using the equation stored in the device; the result of the measurement is shown directly on the display.
 $c \text{ (g/dL)} = 31.4 \times \text{Absorbance (546 nm)}$

Quality Control
 Haemoglobin control material *HEM QS* from Diaglobal GmbH
 Control blood *Para 12 Extend (low, normal, high)* from the Streck company.

Participation in surveys according to the guideline of the German Medical Association for quality assurance of medical laboratory examinations does not apply since the determination is carried out as an individual test in the presence of the patient. Only an internal quality control is required.

Features of analysis performance
Correctness
 The correctness of the analytical results is proven by surveys by the RfB (Reference Institute for Bioanalytics) and control blood measurements.

Method comparison
 A comparison of the Diaglobal test HB 142 (y) with the reference method (x) carried out with human blood showed good conformity.
 Regression analysis
 Passing-Bablok $y = 1.017x - 0.12$
 Lin. Regression $y = 1.015x - 0.112$
 $r = 0.995$
 $n = 37$
 Concentration range: 6.5 - 20.8 g/dL

Precision
 The reproducibility was checked with control blood. Typical VK values are listed in the table below.

In the series N = 20	Average [g/dL]	Standard deviation [g/dL]	VK [%]
Sample 1	5.60	0.07	1.3
Sample 2	11.4	0.10	0.9
Sample 3	16.2	0.15	0.9
From day to day [n = 20]	Average [g/dL]	Standard deviation [g/dL]	VK [%]
Sample 1	5.70	0.09	1.6
Sample 2	11.5	0.13	1.1
Sample 3	16.3	0.20	1.2

Linearity
 The linearity was examined with human and control samples by variation of the sample volume and is given up to a concentration of 30 g/dL.
 The Cusum Test does not show a deviation from the linearity.

Analytical sensitivity
 Lower detection limit: 0.2 g/dL

Analytical specificity
 The determination is highly specific for haemoglobin and its derivatives. COHb and MetHb are also recorded. Severe lipemia and high leukocyte concentrations are disruptive and can simulate excessively high haemoglobin values. Haemolysis does not disturb.

Reference values

	g/dL	mmol/L
Men	14 - 18	8.69 - 11.2
Women	12 - 16	7.45 - 9.93
Pregnant	> 11	> 6.83

Precautionary statements
 Keep away from children
 Cuvettes with cloudy or brownish, discoloured reagent must not be used.

Safety information
 According to Regulation (EEC) Nr.1272/2008 the reagent is not classified as a dangerous mixture. It contains poisonous potassium cyanide in a very small concentration (0.006%), which is under the consideration limit value (0.1%).

A safety data sheet is available on request.
 The reagent does not contain any CMR or endocrine disrupting substances or substances that could cause sensitisation or allergic reactions in the patient or user.
 All serious incidents related to the product must be reported to the manufacturer and the relevant authority.

Disposal instructions
 Do not allow reagent to enter surface water or any sewage system.
 Cuvettes with reagent are classified as hazardous waste. Dispose of it in accordance with local regulations.
 Uncontaminated and completely emptied packaging can be recycled.

Other materials

Item No.	Description
LH 004	Capillaries 10 µL, ring marked
LH 007	Micropipetter
LH 001	Blood lancets
LH 006	Cuvettes stand
LH 011	Alcohol pads for disinfection
LH 009	Cellulose swab