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Reagent for quantitative In-vitro-determination of creatinine in serum / plasma

KRE 121

Order No. KRE 121
Content: 20 tests

Method
Jaffé Reaction, kinetic measurement without deproteinization

Sample material
Serum, plasma or urine
Stability in serum at +2°C to +8°C: 24 hours
Haemolysis disturbs.

Reagents
Content / concentrations:
1. Picric acid (pre-portioned in round cuvettes)
Picric acid 18 mmol/L
2. Buffer solution
Phosphate 10 mmol/L, NaOH 400 mmol/L, pH > 12.0

Safety information
The buffer solution contains 1.6 % sodium hydroxide and is categorised as a dangerous preparation according to EC Directives.
H319: Causes serious eye irritation
H315: Causes skin irritation
Observe the safety advice on the packaging.
A safety data sheet is available on request.¹⁾

Storage and shelf life
Shelf life: The test reagent can be kept at a temperature between +15°C and +25°C until the expiry date indicated on the packaging.

Measurement conditions
Measurement devices: Vario Photometer
Vario Photometer II
Meas. wavelength: 520nm
Temperature: 37°C
2 point measurement Absorption increase during the reaction

Additionally required: Dry block heater 37°C

Measurement range
0.1 – 4.0 mg/dL (8.8 - 354 µmol/L)
In case of exceeding values, dilute the sample 1+5 with physiological saline solution (multiply the result by 6).

Working instructions

- Pipette 1000µL buffer solution into a round cuvette KRE 121. (Save from light. Stability of mixed reagent: 2 hours)
- Incubate this cuvette in the dry block heater for 7 minutes at 37°C.
- Pipette 500µL serum/plasma in this cuvette, mix and incubate in the dry block heater for 1 minute at 37°C.
- Select test <CRE> and insert the round cuvette into the device.
- The result will be displayed in mg/dL after two minutes.

Pipette in single test cuvettes:	
	Analysis
Buffer solution Incubate for 7 minutes at 37°C	1000 µL
Sample	500 µL
Mix and incubate for 1 minute at 37°C	

Calculation
Creatinine concentration:
 $C_{\text{creatinine}} [\text{mg/dL}] = (A2 - A1) \times \text{factor}$

Calculation of creatinine excretion per day:
Creatinine conc. [mg/dL] / 100 g = creatinine / 24 hrs.

Quality assurance
For quality assurance we recommend universal control sera from company Roche, www.roche.de:
PeciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
Order-No.: 05 947 626 190 / 05 947 774 190
Ref.: Roche / Hitachi analyzers,
Method: Jaffé compensated STAT

Reference values²⁾

Serum	mg/dL	µmol/L
Men	below 50 years	0.84 - 1.25
	above 50 years	0.81 - 1.44
Women		0.66 - 1.09

Summary
Serum creatinine is a degradation product which arises due to the decomposition of creatinine phosphate, an important energy storage of the muscle cell.
Patients with a renal disfunction have an elevated serum creatinine level.

Indications / diagnostic significance: ^{2,5)}
Detect a limited glomerular filtration rate (GFR) at
- acute and chronic renal diseases
- pathological urine findings
- hypertonia
- metabolic disorders (diabetes mellitus, hyperuricemia)
- haemodialysis therapy
- extrarenal diseases with diarrhoea ecc.
- screening examinations of symptom-free patients
- monitoring and re-examination of patients with renal diseases.

For the determination of creatinine the economic Jaffé method, which is implemented in diverse variants, has proved best in routine operations. Furthermore, different enzymatic methods exist which show a higher specificity.
The Diaglobal test is based on the Jaffé Reaction. The measurement is performed kinetically.

Measurement principle
Creatinine generates with picric acid in a basic solution a yellow-red colour complex.

Creatinine + Picric acid → Colour complex

The increase of the colour intensity is directly proportional to the creatinine concentration in the sample and is determined photometrically at 520 nm.^{3,4)}

Performance parameters
Specificity / interferences⁶⁾
No interference due to haemoglobin (< 4.0 g/L), bilirubin (< 5 mg/dL), lipaemia (triglycerides < 250 mg/dL), glucose (< 270 mg/dL), urobilinogen (< 40 mg/dL) and ascorbic acid in physiological concentrations (< 25 mg/L). Elevated values because of acetoacetic acid (> 5.0 mmol/L), acetone (> 3.4 mmol/L) and pyruvate (> 1 mmol/L).
Interferences due to pharmaceuticals: falsified elevated values because of cefoxitin (> 70 mmol/L).

Inaccuracy
The reproducibility was checked using human and control samples.

In series [n = 20]	Average [mg/dL]	Standard deviation [mg/dL]	VK [%]
Sample 1 Sample 2	1.14 3.52	0.04 0.09	3.9 2.7
From day to day [n = 20]	Average [mg/dL]	Standard deviation [mg/dL]	VK [%]
Sample 1 Sample 2	1.12 3.50	0.05 0.10	4.4 2.9

Analytic sensitiveness
Lower detection limit: 0.1 mg/dL

Comparison of methods
A comparison of the Diaglobal test KRE 121 (y) and a commercially available test (x) resulted in the following correlation according to the Passing/Bablok ⁷⁾ process:
 $y = 1.038x - 0.141$
 $r = 0.996$
n = 26, concentration range: 0.5 - 4.0 mg/dL

Bibliography

- http://www.diaglobal.de/de/service/downloads/index.html
- Thomas L. Labor und Diagnose. 4th edition. Marburg: Die Medizinische Verlagsgesellschaft 1995: 447
- Bartels H, Böhmer M, Heierli C. Clin Chim Acta 1972; 37:193
- Helger R, Rindfrey H, Hilgenfeldt J. Z Klin Biochem 1974; 12 344
- Sarre H. Nierenkrankheiten. Stuttgart: Georg Thieme, 1959
- Sonntag O. Arzneimittelinterferenzen. Stuttgart, New York: Thieme Verlag, 1985:31
- Passing H, Bablok W. A new biometric procedure for testing the equality of measurements from two different analytical methods. J Clin Chem Clin Biochem. 1983; 21:709-72