**Best. Nr.** MG 013  
**Inhalt:** 2 x 50 mL Farbreagenz  
 Zusätzlich erforderlich:  
 Calcium-Magnesium-Standard Best. Nr. CA ST  
 Diaglobal-Rundküvetten, 40 Stk. Best. Nr. LH 075

**Methode**  
 Xylidylblau<sup>1,2)</sup>

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Contents: 2 x 50 mL Colour reagent

Additionally needed:

Calcium-Magnesium-Standard Order No. CA ST  
Diaglobal-Round Cuvettes, 40 pcs. Order No. LH 075

### Method

Xylydyl blue<sup>1) 2)</sup>

### Sample material

Serum, heparinised plasma or urine

Do not use EDTA-plasma, haemolytical plasma or sera.

Add some drops of HCl to urine, dilute 1+4 with dest. Water.

Stability in serum: at +2°C to +25°C: 7 days

### Reagent

Contents / concentrations of the colour reagent:

Xylydylblue 0.1 mmol/L; GEDTA 40 µmol/L; Potassium carbonate 70 mmol/L; Detergent 1.5%; Tris-buffer 200 mmol/L; Sodium azide < 0.1%; Triton X-100 < 1%, pH 11,0

### Safety information

The reagent contains sodium azide (< 0.1 %) and Triton X-100. Do not swallow and avoid contact with skin and mucous membranes.

If desired a safety data sheet will be provided.<sup>3)</sup>

### Storage and shelf life

The colour reagent can be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging.

### Measurement conditions

Measurement device: VET Photometer Diaglobal

Meas. wavelength: 520 nm

Temperature: Room temperature

### Measurement range

0.10 - 2.10 mmol/L (0.24 - 5.11 mg/dL)

Should values exceed this range, dilute sample 1+1 with physiological saline solution.

Multiply the result by 2.

### Notes

1. Protect reagent from sunlight. Close the bottle after using.
2. The magnesium assay is very sensitive. We recommend the application of disposables.
3. For the "on-site" preparation of plasma we recommend the Diaglobal-Mini centrifuge (Order No. DZ 002)

### Working instructions

Pipette into round cuvettes:			
	Blank	Standard	Analysis
Colour reagent	1500 µL	1500 µL	1500 µL
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix well and measure after 5 minutes at room temperature standard and analysis against blank.

- Select the <MG> test
- Insert cuvette with reagent into photometer (blank value)
- Remove blank cuvette after signal ton
- Insert cuvette with standard into photometer (standard-absorption is shown on the display)
- Remove cuvette
- Insert cuvette with sample (Analysis) into photometer
- Read the result

### Tip

It is possible to prepare several vials for measurement at the same time by pre-pipetting 1500 µL of colour reagent into the appropriate number of Diaglobal round cuvettes.

Shelf life of the pipetted reagents in the vial: 5 hours

### Calculation

Concentration of magnesium in serum/plasma/urine:

$c \text{ [mmol/L]} = 0.8226 \times \text{Abs. analysis} / \text{Abs. standard}$

$c \text{ [mg/dL]} = 2.0 \times \text{Abs. analysis} / \text{Abs. standard}$

Conversion in mg/dL:  $\text{mmol/L} \times 2.4312 = \text{mg/dL}$

### Quality assurance

For quality assurance we recommend our calcium - /magnesium - control solution CA QS.

### Reference values for cattle<sup>4)</sup>

	mmol/L	mg/dL
In serum / plasma	0.8 - 1.3	1.9 - 3.2

### Summary<sup>4)</sup>

Magnesium is a component of many different enzymes which are involved in energy metabolism in the muscle cells. A reduced concentration of magnesium in the blood (hypomagnesaemia) leads to tetanic seizures, while an increase (hypermagnesaemia) leads to flaccid paralysis.

### Evaluation<sup>4)</sup>

Hypermagnesaemia

- Parturient paresis
- Renal insufficiency
- Exsiccosis
- Latrogenic oversupply

Hypomagnesaemia

- Chronic bowel diseases (malabsorption)
- Low-magnesium diet
- Acute pancreatitis
- Cholestasis
- Renal insufficiency
- Hypercalcaemia
- Hypocalcaemic ketoacidosis
- Hyperthyroidism
- Tetany (cow, sheep)

### Measurement principle<sup>1) 2)</sup>

In alkaline solution with xylydyl blue, magnesium forms a colour complex which is measured at 520 nm. Calcium ions are masked with Glycol Ethylenediaminetetraacetic Acid (GEDTA).

The assessment is carried out using a standard.

$\text{Mg}^{2+} + \text{Xylydyl blue} \rightarrow \text{Magnesium colour complex}$

### Performance parameters

#### Specificity / interferences

No interference caused by calcium ions.

Bilirubin up to 20 mg/dL and lipidaemia do not interfere.

Do not use haemolytic serums as erythrocytes have high magnesium concentrations.

#### Inaccuracy

The reproducibility was checked using human and control samples.

In series [n = 20]	Average [mmol/L]	Standard deviation [mmol/L]	VK [%]
Sample 1	0.88	0.02	2.3
Sample 2	1.15	0.02	2.1
Sample 3	1.45	0.02	1.6
From day to day [n = 20]	Average [mmol/L]	Standard-deviation [mmol/L]	VK [%]
Sample 1	0.67	0.04	5.8
Sample 2	0.93	0.05	5.4
Sample 3	1.38	0.05	3.9

### Analytic sensitiveness

Lower detection limit: 0.08 mmol/L (0.2 mg/dL)

### Comparison of methods

A comparison of the Diaglobal test MG 013 (y) and a commercially available test (x) resulted in the following correlation according to the linear regression:

$$y = 0,982x + 0,02$$

$$r = 0,993$$

n = 26, concentration range: 0.5 - 2.0 mmol/L

### Information on disposal

Waste code number 180106:

Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system.

Dispose of in accordance with official regulations.

Non-contaminated and completely empty packaging can be recycled.

Non-contaminated and completely empty packaging can be recycled.

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