

Survey: KS4/19
month: May 2019
Participant-No: 0004965
date of issue: 18.05.19



Referenzinstitut für Bioanalytik

diaglobal GmbH
Dr. Rainer Schlegel
Haus 41
Köpenicker Str. 325
12555 Berlin



Survey directors

Prof. Dr. C. Knabbe
Prof. Dr. Dr. K.P. Kohse
Prof. Dr. M. Neumaier

Head of RfB

Dr. W. J. Geilenkeuser
Dr. A. Kessler

Bonn, 14. June 2019

Certificate

We confirm that you have participated in the survey for clinical chemical analytes in serum.

You have met the requirements of the survey for the following analytes:

bilirubin, total	(2)	cholesterol	(1)	glucose	(3)
urea	(1)	creatinine	(1)	lactate	(2)
triglycerides	(2)	calcium	(3)	magnesium	(4)
CK	(1)	AST	(1)	ALT	(1)

This certificate is valid until the end of November 2019.

The number in parentheses characterizes the analytical method used.
The assignment of the number to the respective method and/or the respective instrument is to be taken from the total evaluation.

Survey: **KS4/19**
 month: **May 2019**
 Participant-No: **0004965**



Listing and Evaluation of all your results

Explanations

Certificate

A certificate is issued (given) for an analyte only if the basis for an evaluation of the accuracy is given by the guidelines of the German Medical Association and/or if an evaluation is possible in analogy to these guidelines (see comments on the evaluation) and, if both results for an analyte are within the given acceptance limits. (marked as '+' below C)

Certificate of participation

In the participation certificate all analytes which are included in the list on this page are listed.
 If all analytes are listed on the certificate no participation certificate is printed.

Legend:	C=Certification, M=No of method, R=your result, D=difference (R-T) Dmax= maximum allowable amount of difference in measurement, partly given by annex 1 of the guideline of the BAK (Dt. Arzteblatt 111, Heft 38, 19.9.2014). T = target value, either reference method value or assigned value, LL UL = lower resp. upper limit	Certification:	+ = fulfilled (quotient D/Dmax <= 1.0) - = not fulfilled (quotient D/Dmax > 1.0 ± = certification cancelled because of technical and/or analytical reasons
---------	--	----------------	--

	C	M	R	D/Dmax	T	LL	UL							
bilirubin, total [mg/dl]	+	2	A	3.16	0.11	3.08	2.40	3.76						
			B	1.79	-0.03	1.80	1.40	2.20						
cholesterol [mg/dl]	+	1	A	113	0.12	111	96.7	126						
			B	133	0	133	115	151						
glucose [mg/dl]	+	3	A	92.4	-0.13	94.2	80.0	109						
			B	117	-0.06	118	100	136						
urea [mg/dl]	+	1	A	66.8	0.50	60.7	48.5	72.8						
			B	119	0.08	117	93.6	141						
creatinine [mg/dl]	+	1	A	5.30	-0.16	5.48	4.38	6.58						
			B	2.06	0.52	1.86	1.48	2.24						
lactate [mmol/l]	+	2	A	3.54	0.30	3.35	2.74	3.97						
			B	4.18	0.39	3.90	3.18	4.61						
triglycerides [mg/dl]	+	2	A	228	0.34	216	181	251						
			B	126	0.25	121	101	141						
calcium [mmol/l]	+	3	A	2.47	0.26	2.41	2.16	2.65						
			B	1.95	0.05	1.94	1.74	2.14						
magnesium [mmol/l]	+	4	A	1.41	-0.19	1.45	1.23	1.67						
			B	0.960	0.00	0.959	0.815	1.11						
CK [U/l]	+	1	A	530	-0.23	556	444	667						
			B	358	0.32	336	268	404						
AST [U/l]	+	1	A	252	0.13	245	193	297						
			B	345	0.02	343	270	416						
ALT [U/l]	+	1	A	186	-0.52	209	165	254						
			B	181	-0.10	185	146	224						